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English

Thank you for your purchase of an FH-G1 MEDICAL SLIDE HOLDER.

The FH-G1 is designed specifically for scanning slides glass (for microscope) in Nikon's SUPER COOLSCAN 4000 ED (4000 ED), COOLSCAN IV ED, LS-2000, COOLSCAN III to correspond to the slide holder.

Slide glass:
26mm (width) × 76mm (length) × 0.8 ~ 1.5mm (thickness)

Cover glass:
25mm below (width) × 60mm below (length) × 0.18mm below (thickness)

Slide glass and cover glass are sealed with adhesive.

Thickness (including adhesive):
2mm or less (includes increase in thickness due to adhesive)

Do not use opaque samples.

Before Scanning

1. Hold the FH-G1 with the Nikon logo facing up. Gently slip the prepared slide, with the cover uppermost, into the slide insertion slot (Figure B).

2. Adjust the position of the prepared slide until the portion you wish to scan is within the scan area (⑤ in Figure A, indicated by the scanning guide marks ④).

Caution

- After positioning the prepared slide, handle the holder carefully to prevent the slide from slipping. Although the slide is held in place by the retaining clips inside the slot, jolting the holder could cause the slide to slip out of position.

3. Insert the MA-20 or MA-20 (S) SLIDE MOUNT ADAPTER into the scanner and check that both the scanner and computer are on.

4. Holding the FH-G1 with the Nikon logo facing up, slide the holder into the MA-20 or MA-20 (S) film slot in the direction indicated by the arrow (Figure C).

Notices

- The information contained in this manual is subject to change without notice.

- Nikon does not accept liability for damage to slides glass (for microscope) resulting from the use of this product.

For Your Safety

Be sure to read this manual thoroughly before using the product. After reading, store this manual where it will be read by all who use the product.

Precautions for Use

- If the product is dropped or subject to other forms of physical shock, the product or prepared slide could break, resulting in injury.
- Prepared slides of a thickness greater than the limit imposed by slide insertion guides at the holder entrance can not be used. Attempting to force a prepared slide of greater thickness into the holder could cause the holder or slide to break, resulting in injury.

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After Scanning

After completing the scan, reverse the order of the steps listed above under "Before Scanning" to remove the holder from the scanner and take the prepared slide from the holder.

Caution

- You may find that the slide sticks slightly when removed from the holder. This is normal, and the slide can be removed from the holder by slipping it out slowly.

Using the FH-G1

Slides Glass (for microscope)

This product can be used with slides of the dimensions given below.

Slide: 26mm (W) × 76mm (L), 0.8–1.5mm in thickness

Cover: 25mm (W) × 60mm (L) or less, thickness of 0.18mm or less

Total thickness (after slide and cover are sealed with cement):

2mm or less (includes increase in thickness due to cement)

Do not use opaque samples.

Maintenance

Should the product become dirty, clean it gently with a soft cloth.

Do not expose the product to volatile chemicals such as alcohol, thinner, or benzene, as this could cause discoloration or other deterioration.

Specifications

Effective scan area:

22.9 mm × 35 mm (area indicated by scanning guide marks)

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